



## Hornet IPG

### Minor Mechanical Shock DVT Protocol

Effective Date:

Rev:

Pg. 1 of 10

#### 1 Purpose

This protocol prescribes methods and records results necessary to verify minor mechanical shock testing, minor mechanical shock testing, and drop and impact testing meets requirements for the 3025 Hornet IPG.

#### 2 Scope

This document references specific standards and device requirements to be tested by an external laboratory. This document provides records to ensure testing is performed to the required standards, and results are reviewed at the completion of testing

#### 3 References

Document No.	Title
ISO 14708-3: 23.2	Implants for surgery – Active implantable medical devices – Part 1: General requirements for safety, marking and for information to be provided by the manufacturer, Part 23.2, Protection of the ACTIVE IMPLANTABLE MEDICAL DEVICE from mechanical forces
EN 60068-2-47	Environmental testing – Part 2-47: Tests – Mounting of specimens for vibration, impact, and similar dynamic tests
EN 60068-2-64	Environmental testing – Part 2-64 – Test Fh: Vibration, broadband random and guidance
ISO 14708-3: 23.7	Implants for surgery – Active implantable medical devices – Part 1: General requirements for safety, marking and for information to be provided by the manufacturer, Part 23.7, Protection of the ACTIVE IMPLANTABLE MEDICAL DEVICE from mechanical forces
EN 60068-2-27	Environmental testing – Part 2-27: Tests – Test Ea and guidance: Shock
ASTM D5276	Standard Test Method for Drop Test of Loaded Containers by Free Fall
ASTM D4332	Standard Practice for Conditioning Containers, Packages, or Packaging Components for Testing

#### 4 Appendices

Appendix:	Title
A	Minor Mechanical Shock Testing Request Form
B	Functional Test Record
C	Additional Notes Area (if required)

CONFIDENTIAL

All information contained in this document is confidential and is the sole property of Med-Ally Inc. Any reproduction in part or whole without the written permission of Med-Ally Inc. is prohibited



Hornet IPG		
Minor Mechanical Shock DVT Protocol		
Effective Date:	Rev:	Pg. 2 of 10

5 Definitions

Abbreviation or Term	Definition
DVT	Design Verification Test
IPG	Implantable Pulse Generator



Hornet IPG		
Minor Mechanical Shock DVT Protocol		
Effective Date:	Rev:	Pg. 3 of 10

**6 Testing Protocol**

- 6.1 Devices will be tested to the following requirements:
  - 6.1.1 Minor Mechanical Shock Testing:
    - 6.1.1.1 Shock Shape: Half sine or Haversine
    - 6.1.1.2 Severity: Peak Acceleration: 5000m/s^2 (500g)
    - 6.1.1.3 Duration of Shock: 1ms
    - 6.1.1.4 Shock Details: Six total; One shock in each direction along each axis XYZ
    - 6.1.1.5 Functional Testing: Device meets all functional testing requirements after completion of Minor Mechanical Shock testing

6.2 Testing Device Drawing Number: \_\_\_\_\_

6.3 Functional Test Work Instruction: \_\_\_\_\_

**6.4 Information for Outsourcing:**

- 6.4.1 Complete Appendix A, Minor Mechanical Shock Testing Request Form. Submit the form to the anticipated vendor to obtain quoting information and attach the provided quote to this protocol.
- 6.4.2 Print a copy of the applicable assembly drawing(s) to include with the request form
- 6.4.3 Print a photograph of the completed devices to include with the request form

**Signature:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**6.5 Approval:**

A representative from QA must review and approve the specification information and submission for prior to shipment.  
Quality Approval to Execute Testing:

**Signature:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**6.6 Sample Generation**

- 6.6.1 Describe the origin of the samples used. If new samples where created for this test, describe any deviations, if applicable, that may impact testing:

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

CONFIDENTIAL

All information contained in this document is confidential and is the sole property of Med-Ally Inc. Any reproduction in part or whole without the written permission of Med-Ally Inc. is prohibited



Hornet IPG		
Minor Mechanical Shock DVT Protocol		
Effective Date:	Rev:	Pg. 4 of 10


Initial: \_\_\_\_\_ Date: \_\_\_\_\_

Record Device Serial Numbers Below:




## Hornet IPG

### Minor Mechanical Shock DVT Protocol

Effective Date:

Rev:

Pg. 5 of 10

---

---

---

---

---

Initial: \_\_\_\_\_ Date: \_\_\_\_\_

#### 6.7 Lab Testing:

**NOTE: Any testing laboratory that can meet the testing requirements may be used.**

- 6.7.1 Complete a pre-testing functional test prior to shipment. Record the results in Appendix B: Functional Test Record

Initial: \_\_\_\_\_ Date: \_\_\_\_\_

- 6.7.2 Ensure the device battery is less than 30% charged prior to shipping.

Initial: \_\_\_\_\_ Date: \_\_\_\_\_

- 6.7.3 Package the devices for shipment as to avoid damage in transit. Include the request form, device drawing and photograph, and any other relevant work order information with the devices.

Initial: \_\_\_\_\_ Date: \_\_\_\_\_

- 6.7.4 Prior to shipping review all paperwork for completeness. Include a review by a quality representative for verification.

Signature: \_\_\_\_\_ Date: \_\_\_\_\_

(QA) Signature: \_\_\_\_\_ Date: \_\_\_\_\_

- 6.7.5 Upon Receipt:

- 6.7.5.1 Review the provided reports for completeness:

Initial: \_\_\_\_\_ Date: \_\_\_\_\_

#### CONFIDENTIAL

All information contained in this document is confidential and is the sole property of Med-Ally Inc. Any reproduction in part or whole without the written permission of Med-Ally Inc. is prohibited



Hornet IPG		
Minor Mechanical Shock DVT Protocol		
Effective Date:	Rev:	Pg. 6 of 10

6.7.5.2 Complete a Functional Test, and record the results for each device on Appendix B, Functional Test Record. Attach copies of any functional test reports generated by testing.

Initial: \_\_\_\_\_ Date: \_\_\_\_\_

6.7.5.3 Attach all lab data to this protocol:

Initial: \_\_\_\_\_ Date: \_\_\_\_\_  
Notes:

---

---

---

---

---



Hornet IPG		
Minor Mechanical Shock DVT Protocol		
Effective Date:	Rev:	Pg. 7 of 10

7 Approvals

7.1 Verify testing results meet drawing requirements

Initial\_\_\_\_\_ Date\_\_\_\_\_

7.2 Quality Approval

- 7.2.1 Review Protocol
- 7.2.2 Review Lab Testing Records
- 7.2.3 Review Appendix B: Functional Test Record
- 7.2.4 Review Appendix C: Additional Notes Area (if applicable)
- 7.2.5 Ensure Testing Results are attached

Signature:\_\_\_\_\_ Date:\_\_\_\_\_

7.3 Other Approval (If required):

Signature:\_\_\_\_\_ Date:\_\_\_\_\_

7.4 Notes (if required):

---

---

---

---

---

---

---

---

---

---



Hornet IPG

Minor Mechanical Shock DVT Protocol

Effective Date:

Rev:

Pg. 8 of 10

## Appendix A: Minor Mechanical Shock Testing Request Form

Point of Contact: \_\_\_\_\_ Ph. \_\_\_\_\_

**Bill To:** Med-Ally  
2040 Bushy Park Rd.  
CIMC N. Bldg 6  
Goose Creek, SC 29445

**Ship to:** Med-Ally  
2040 Bushy Park Rd.  
CIMC N. Bldg 6  
Goose Creek, SC 29445

**Project Name:** Hornet IPG

**Product Weight:** \_\_\_\_\_ oz each; \_\_\_\_\_ oz total shipped weight

Quantity to test: \_\_\_\_\_

**Description of Contents:** Medical device; Internal Pulse generator consisting of active electronics with Lithium Ion batteries contained in a welded titanium enclosure, with external connections accessible in the attached epoxy header.

**Special Handling:** BATTERY CHEMISTRY, WEIGHT battery in each device; No other drugs, biologics or hazardous materials included.

**Battery Considerations:** NOTE: Westpak requires batteries to be UN 38.3 certified with a state of charge <30%, and include MSDS

**Max Temperature Exposure:** \_\_\_\_\_ **Min Temp Exposure:** \_\_\_\_\_

**Testing Required:** Vibration Mechanical Force Testing in accordance with ISO 14708-3: 23.7 and EN 60068-2-27, reference method below:

- a) Shock Shape: Half sine or Haversine
- b) Severity: Peak Acceleration: 5000m/s<sup>2</sup> (500g)
- c) Duration of Shock: 1ms
- d) Shock Details: Six Total: One shock in each direction along each axis XYZ
- e) Return devices to Med-Ally upon test completion for Functional Verification
- f) Include verification of test methods, dates, times and associated information with returned parts.

QA Signature: \_\_\_\_\_ Date: \_\_\_\_\_

CONFIDENTIAL

All information contained in this document is confidential and is the sole property of Med-Ally Inc. Any reproduction in part or whole without the written permission of Med-Ally Inc. is prohibited





Hornet IPG		
Minor Mechanical Shock DVT Protocol		
Effective Date:	Rev:	Pg. 9 of 10

**Appendix B: Functional Test Record**

Circle the functional test being completed below:

Pre-Mechanical Shock Testing

Post-Mechanical Shock Testing

Serial #	Pass/Fail	Notes:

Signature: \_\_\_\_\_ Date: \_\_\_\_\_ Page \_\_\_\_ of \_\_\_\_



Hornet IPG		
Minor Mechanical Shock DVT Protocol		
Effective Date:	Rev:	Pg. 10 of 10

**Appendix C: Additional Notes (if required). Notes may be typed or hand written:**

Signature: \_\_\_\_\_ Date: \_\_\_\_\_ Page \_\_\_\_ of \_\_\_\_